

Straumann® BLX Implant System

INTERNATIONAL EDITION

Product Catalog 2022/2023

Roxolid® SLA®

Product	Platform	Image	Material	Length	Art.-No.
Straumann® BLX Implant Ø 3.5 mm	RB – Regular Base Ø 3.5 mm		Roxolid®	8 mm	061.3508
				10 mm	061.3510
				12 mm	061.3512
				14 mm	061.3514
				16 mm	061.3516
				18 mm	061.3518
Straumann® BLX Implant Ø 3.75 mm	RB – Regular Base Ø 3.75 mm			6 mm	061.4506
				8 mm	061.4508
				10 mm	061.4510
				12 mm	061.4512
				14 mm	061.4514
				16 mm	061.4516
Straumann® BLX Implant Ø 4.0 mm	RB – Regular Base Ø 4.0 mm			18 mm	061.4518
				6 mm	061.5506
				8 mm	061.5508
				10 mm	061.5510
				12 mm	061.5512
				14 mm	061.5514
Straumann® BLX Implant Ø 4.5 mm	RB – Regular Base Ø 4.5 mm			16 mm	061.5516
				18 mm	061.5518
				6 mm	061.6506
				8 mm	061.6508
				10 mm	061.6510
				12 mm	061.6512
				14 mm	061.6514
				16 mm	061.6516
				18 mm	061.6518

Surgery

RB Closure Cap

Product	Image	Description	Height	Material	Art.-No.
Closure Caps					
RB Closure Cap			0,4 mm	Ti	064.41005

Roxolid® SLA®

Product	Platform	Image	Material	Length	Art.-No.
Straumann® BLX Implant Ø 5.0 mm	WB – Wide Base Ø 5.0 mm		Roxolid®	6 mm	061.7506
				8 mm	061.7508
				10 mm	061.7510
				12 mm	061.7512
				14 mm	061.7514
				16 mm	061.7516
				18 mm	061.7518
Straumann® BLX Implant Ø 5.5 mm	WB – Wide Base Ø 5.5 mm			6 mm	061.8506
				8 mm	061.8508
				10 mm	061.8510
				12 mm	061.8512
				14 mm	061.8514
				16 mm	061.8516
Straumann® BLX Implant Ø 6.5 mm	WB – Wide Base Ø 6.5 mm			6 mm	061.9506
				8 mm	061.9508
				10 mm	061.9510
				12 mm	061.9512
				14 mm	061.9514
				16 mm	061.9516

Surgery


WB Closure Cap

Product	Image	Description	Height	Material	Art.-No.
Closure Caps					
WB Closure Cap			0,5 mm	Ti	064.8102S

Healing Abutments

Product	Image	Description	Material	Art.-No.
Healing Abutments, Crown				
RB/WB Healing Abutment, for final abutments Ø 3.8 mm		Crown, Ø 4 mm, gingiva height 1.5 mm, abutment height 2 mm, (total 3.5 mm)	Ti	064.4202S
		Crown, Ø 4 mm, gingiva height 1.5 mm, abutment height 4 mm, (total 5.5 mm)		064.4203S
		Crown, Ø 4 mm, gingiva height 2.5 mm, abutment height 2 mm, (total 4.5 mm)		064.4204S
		Crown, Ø 4 mm, gingiva height 2.5 mm, abutment height 4 mm, (total 6.5 mm)		064.4205S
		Crown, Ø 4 mm, gingiva height 3.5 mm, abutment height 2 mm, (total 5.5 mm)		064.4206S
		Crown, Ø 4 mm, gingiva height 3.5 mm, abutment height 4 mm, (total 7.5 mm)		064.4207S
RB/WB Healing Abutment, for final abutments Ø 4.5 mm		Crown, Ø 5 mm, gingiva height 1.5 mm, abutment height 2 mm, (total 3.5 mm)	Ti	064.4212S
		Crown, Ø 5 mm, gingiva height 1.5 mm, abutment height 4 mm, (total 5.5 mm)		064.4213S
		Crown, Ø 5 mm, gingiva height 2.5 mm, abutment height 2 mm, (total 4.5 mm)		064.4214S
		Crown, Ø 5 mm, gingiva height 2.5 mm, abutment height 4 mm, (total 6.5 mm)		064.4215S
		Crown, Ø 5 mm, gingiva height 3.5 mm, abutment height 2 mm, (total 5.5 mm)		064.4216S
		Crown, Ø 5 mm, gingiva height 3.5 mm, abutment height 4 mm, (total 7.5 mm)		064.4217S

Healing Abutment Dimensions



Healing Abutments

Product	Image	Description	Material	Art.-No.
RB/WB Healing Abutment, for final abutments Ø 6.0 mm		Crown, Ø 6.5 mm, gingiva height 1.5 mm, abutment height 2 mm, (total 3.5 mm)	Ti	064.4222S
		Crown, Ø 6.5 mm, gingiva height 1.5 mm, abutment height 4 mm, (total 5.5 mm)		064.4223S
		Crown, Ø 6.5 mm, gingiva height 2.5 mm, abutment height 2 mm, (total 4.5 mm)		064.4224S
		Crown, Ø 6.5 mm, gingiva height 2.5 mm, abutment height 4 mm, (total 6.5 mm)		064.4225S
		Crown, Ø 6.5 mm, gingiva height 3.5 mm, abutment height 2 mm, (total 5.5 mm)		064.4226S
		Crown, Ø 6.5 mm, gingiva height 3.5 mm, abutment height 4 mm, (total 7.5 mm)		064.4227S
		Crown, Ø 7.5 mm, gingiva height 2.5 mm, abutment height 2 mm, (total 4.5 mm)		064.4234S
		Crown, Ø 7.5 mm, gingiva height 2.5 mm, abutment height 4 mm, (total 6.5 mm)		064.4235S
		Crown, Ø 7.5 mm, gingiva height 3.5 mm, abutment height 2 mm, (total 5.5 mm)		064.4236S
		Crown, Ø 7.5 mm, gingiva height 3.5 mm, abutment height 4 mm, (total 7.5 mm)		064.4237S
Healing Abutments, Bridge				
RB/WB Healing Abutment, for final abutments Ø 4.5 mm		Bridge/Bar, Ø 5 mm, gingiva height 1.5 mm, abutment height 2 mm, (total 3.5 mm)	Ti	064.4232S
		Bridge/Bar, Ø 5 mm, gingiva height 1.5 mm, abutment height 4 mm, (total 5.5 mm)		064.4233S

Always communicate the selected healing abutment to the dental lab

The BLX Healing Abutment selection determines the available options for the final restoration. Therefore, always communicate the implant diameter and the selected healing abutment to the restoring dentist or dental lab in order to facilitate the appropriate final abutment selection.

For more information refer to *Straumann® BLX Implant System – Healing Abutment selection, Quick Guide* (704658/en).

Straumann® BLX Implant System – Healing Abutment selection Quick Guide

The total height of a BLX Healing Abutment consists of two parts:

Gingiva Height (GH): Emerging part of the healing abutment, available in heights of 0.5 mm to 3.5 mm. Uniform profile throughout the BLX prosthetic portfolio.

Abutment Height (AH): Cylindrical part of the healing abutment, available in two heights: 2 mm (low) or 4 mm (high).

Always communicate the implant diameter and the selected healing abutment to the dental lab!

Please consider, that the BLX Healing Abutment selection determines the available options for the final restoration.

Slim profile single crown

RB/WB Healing Abutment



After healing, no wide profile crown (WB) nor bridge restoration possible, even on the back side.

Note: Use Healing Abutment Ø5.0, GH 1.5 mm for Variobase® AL and RB/ WB CARISTY® Abutments.

Wide profile single crown, molar

WB Healing Abutment



Note: For WB BLX implants.

Bridge

RB/WB Healing Abutment for Bridge/Bar



Note: Use this healing abutment if it has not yet been determined whether the restoration will consist of single crowns or a bridge.

For more information refer to the *Straumann® BLX Implant System Consistent Emergency Pro®/im® and Healing Abutment selection* (704658/en).

straumann