



Tapered Pro Conical

Catalog & surgical manual

Tapered Conical


Implant portfolio

Tapered Pro Conical

						
Body diameter	3.3mm	3.8mm	3.8mm	4.2mm	4.6mm	5.2mm
Prosthetic connection	 Narrow	 Narrow	 Regular	 Regular	 Regular	 Regular
Platform switch	0.3mm	0.3mm	0.4mm	0.4mm	0.6mm	0.9mm
Laser-Lok zone	1.8mm	1.8mm	1.8mm	1.8mm	1.8mm	1.8mm
Apical diameter	2.4mm	2.8mm	2.8mm	3.0mm	3.2mm	3.3mm
9.0mm length	CTA3309	CTA3809	CTA3809R	CTA4209	CTA4609	CTA5209
10.5mm length	CTA3310	CTA3810	CTA3810R	CTA4210	CTA4610	CTA5210
12.0mm length	CTA3312	CTA3812	CTA3812R	CTA4212	CTA4612	CTA5212
15.0mm length	CTA3315	CTA3815	CTA3815R	CTA4215	CTA4615	CTA5215
18.0mm length	CTA3318	CTA3818	CTA3818R	CTA4218	CTA4618	—

Tapered Short Conical

			
Body diameter	4.2mm	4.6mm	5.2mm
Prosthetic connection	 Regular	 Regular	 Regular
Platform switch	0.6mm	0.8mm	1.0mm
Laser-Lok zone	1.8mm	1.8mm	1.8mm
Apical diameter	3.3mm	3.7mm	4.2mm
7.5mm length	CTA4207	CTA4607	CTA5207



Healing & temporary

Cover screws, healing caps & temporary abutments

Conical Cover Screw



Part number	Diameter
CNCC	Narrow
CRCC	Regular

Regular Healing Caps



Part number	Prosthetic connection	Gingival height	Gingival diameter
CNHA2	Narrow	2.0mm	3.0mm
CNHA4		4.0mm	3.0mm
CNHA6		6.0mm	3.0mm
CRHA2	Regular	2.0mm	3.8mm
CRHA4		4.0mm	3.8mm
CRHA6		6.0mm	3.8mm

Wide Healing Caps



CNHA4W	Narrow	4.0mm	4.8mm
CRHA4W	Regular	4.0mm	5.3mm
CRHA6W		6.0mm	5.3mm

Extra Wide Healing Caps



CRHA4EW	Regular	4.0mm	5.8mm
CRHA6EW		6.0mm	5.8mm

Hand-tighten (10-15 Ncm) with the .050" (1.25mm) Hex Driver. Titanium alloy.

Conical Temporary Abutment, Engaging



Part number	Platform	Gingival height
CNTEA	Narrow	1.5mm
CNTEA3		3.0mm
CNTEA4		4.0mm
CRTEA	Regular	1.5mm
CRTEA3		3.0mm
CRTEA4		4.0mm

Preparable. Includes abutment screw. Titanium alloy. Tighten to 20Ncm.

Conical Temporary Abutment, Non-engaging



Part number	Platform	Gingival height
CNTNA	Narrow	1.5mm
CNTNA3		3.0mm
CNTNA4		4.0mm
CRTNA	Regular	1.5mm
CRTNA3		3.0mm
CRTNA4		4.0mm

Preparable. Includes abutment screw. Titanium alloy. Tighten to 20Ncm.